

# Animal Water Feeding Bottles

## FOR FEEDING WATER TO EXPERIMENTAL ANIMALS

Orchid Scientific produces a wide range of water bottles. The bottles are made of non-transparent Polypropylene (PP), transparent Polycarbonate (PC) & Polysulfone (PS) with Silicon Bush OR SS Cap.



**Designed To Feed  
Water To Experimental  
Animals**



## FEATURES:




- Autoclavable\*
- Leakage proof- No water leakage from bottle
- Assured water delivery- Each bottle cap is assured for proper water delivery
- Consistent thickness
- Wide mouth to aid filling and emptying, whilst also ensuring effective washing
- Volume graduations to aid volume measurements
- Vacuum style, without ball

**Options of Cap:** Silicon bush with Press fit Style OR S.S. press fit Cap with silicon gasket Press fit Style

## MICE WATER BOTTLES

(Polypropylene, Polycarbonate, Polysulfone)

(Vacuum style, without ball with Silicon Bush / SS Cap Style)

Image	Model No	Capacity	Bottle Type Vacuum Style Without Ball	Cap	Sipper tube type
<b>Polypropylene Bottle for Mice</b>					
 <p>POLYPROPYLENE 150 ML ROUND SHAPE WITH SILICON BUSH MODEL NO: SI-15P</p>	SI-15P	150ml Round Shape	Polypropylen	Silicon bush with S.S. Cap Press fit Style	S.S. 304 Nozzle
 <p>POLYPROPYLENE 150 ML ROUND SHAPE WITH SS CAP MODEL NO: SS-15P</p>	SS-15P	150ml Round Shape	Polypropylene	S.S. press fit Cap with silicon gasket Press fit Style	S.S. 304 Nozzle
<b>Polycarbonate Bottle for Mice</b>					
 <p>POLYCARBONATE 300 ML SQUARE SHAPE WITH SS CAP MODEL NO: SS-30C</p>	SS-30C	300ml Round Shape	Polycarbonate	S.S. press fit Cap with silicon gasket Press fit Style	S.S. 304 Nozzle

Note: Above dimensions are in mm

Image	Model No	Capacity	Bottle Type Vacuum Style Without Ball	Cap	Sipper tube type
<b>Polysulfone Bottle for Mice</b>					
<p>POLYSULFONE 300 ML SQUARE SHAPE WITH SS CAP MODEL NO: SS-30S</p>	SS-30S	300ml Square Shape	Polysulfone	S.S. Cap with silicon gasket Press fit Style	S.S. 304 Nozzle



## RAT WATER BOTTLES

(Polypropylene, Polycarbonate, Polysulfone)



(Vacuum style, without ball with Silicon Bush / SS Cap Style)

Image	Model No	Capacity	Bottle Type Vacuum Style Without Ball	Cap	Sipper tube type
<b>Polypropylene Bottle for Rat</b>					
<p>POLYPROPYLENE 300 ML ROUND SHAPE WITH SILICON BUSH MODEL NO: SI-30P</p>	SI-30P	300ml Round Shape	Polypropylene	Silicon bush with S.S. Cap Press fit Style	S.S. 304 Nozzle

Note: Above dimensions are in mm

Image	Model No	Capacity	Bottle Type Vacuum Style Without Ball	Cap	Sipper tube type
<b>Polypropylene Bottle for Rat</b>					
 <p data-bbox="165 891 504 981">POLYPROPYLENE 300 ML ROUND SHAPE WITH SS CAP MODEL NO: SS-30P</p>	SS-30P	300ml Round Shape	Polypropylene	S.S. press fit Cap with silicon gasket Press fit Style	S.S. 304 Nozzle
<b>Polycarbonate Bottle for Rat</b>					
 <p data-bbox="181 1563 504 1637">POLYCARBONATE 300 ML SQUARE SHAPE WITH SS CAP MODEL NO: SS-30C</p>	SS-30C	300ml Square Shape	Polycarbonate	S.S. press fit Cap with silicon gasket Press fit Style	S.S. 304 Nozzle

Note: Above dimensions are in mm

Image	Model No	Capacity	Bottle Type Vacuum Style Without Ball	Cap	Sipper tube type
<b>Polysulfone bottle for Rat</b>					
 <p data-bbox="212 728 483 795">POLYSULFONE 300 ML SQUARE SHAPE WITH SS CAP MODEL NO: SS-30S</p>	SS-30S	300ml Square Shape	Polysulfone	S.S. Cap with silicon gasket Press fit Style	S.S. 304 Nozzle
 <p data-bbox="212 1191 483 1258">POLYSULFONE 500 ML SQUARE SHAPE WITH SS CAP MODEL NO: SS-50S</p>	SS-50S	500ml Square Shape	Polysulfone	S.S. Cap with silicon gasket Press fit Style	S.S. 304 Nozzle

Note: Above dimensions are in mm

## RABBIT & GUINEA PIG WATER BOTTLES

(Polypropylene, Polycarbonate, Polysulfone)

(Silicon Bush - Vacuum Style Without Ball OR Ball Type / SS Cap style - Vacuum Style Without Ball)

Image	Model No	Capacity	Bottle Type Vacuum Style Without Ball	Cap	Sipper tube type
<b>Polypropylene Bottle for Rabbit / Guinea Pig</b>					
<p>POLYPROPYLENE 650 ML ROUND SHAPE WITH SILICON BUSH SI65P-ST</p>	SI65P-ST (Vacuum Type)	650 ml Round Shape	Polypropylene Vacuum Style Without Ball	Silicon bush with S.S. Cap Press fit Style	Straight Type S.S. 304 Nozzle
<p>POLYPROPYLENE 650 ML ROUND SHAPE WITH SILICON BUSH SI65P-ST</p>	SI65P-ST (Ball Type)	650 ml Round Shape	Polypropylene Ball Type	Silicon bush with S.S. Cap Press fit Style	Straight Type S.S. 304 Nozzle
<p>POLYPROPYLENE 650 ML ROUND SHAPE WITH SILICON BUSH SI65P-BT</p>	SI65P-BT (Vacuum Type)	650 ml Round Shape	Polypropylene Vacuum Style Without Ball	Silicon bush with S.S. Cap Press fit Style	Bent Type S.S. 304 Nozzle

Note: Above dimensions are in mm

Image	Model No	Capacity	Bottle Type Vacuum Style Without Ball	Cap	Sipper tube type
<b>Polypropylene Bottle for Rabbit / Guinea Pig</b>					
<p>POLYPROPYLENE 650 ML ROUND SHAPE WITH SILICON BUSH SI65P-BT</p>	SI65P- BT (Vacuum Type)	650 ml Round Shape	Polypropylene Vacuum Style Without Ball	Silicon bush with S.S. Cap Press fit Style	Bent Type S.S. 304 Nozzle
<p>POLYPROPYLENE 650 ML ROUND SHAPE WITH SILICON BUSH SI65P-BT</p>	SI65P-BT (Ball Type)	650 ml Round Shape	Polypropylene Ball Type	Silicon bush with S.S. Cap Press fit Style	Bent Type S.S. 304 Nozzle
<p>POLYPROPYLENE 650 ML ROUND SHAPE WITH SILICON BUSH SI65P-BT</p>	SI65P-BT (Vacuum Type)	650 ml Round Shape	Polypropylene Vacuum Style Without Ball	Silicon bush with S.S. Cap Press fit Style	Bent Type S.S. 304 Nozzle

Note: Above dimensions are in mm




Image	Model No	Capacity	Bottle Type Vacuum Style Without Ball	Cap	Sipper tube type
<b>Polypropylene Bottle for Rabbit / Guinea Pig</b>					
 <p data-bbox="225 734 469 801">POLYPROPYLENE 650 ML ROUND SHAPE WITH SS CAP SS65P-BT</p>	SS65P-BT	650ml Round Shape	Polypropylene Vacuum Style Without Ball	S.S. Cap with silicon gasket Press fit Style	Bent Type S.S. 304 Nozzle
<b>Polycarbonate Bottle for Rabbit / Guinea Pig</b>					
 <p data-bbox="236 1308 464 1364">POLY CARBONATE 500 ML SQUARE SHAPE WITH SS CAP SS-50C-ST</p>	SS50C-ST	500 ml Square Shape	Polycarbonate Vacuum Style Without Ball	S.S. Cap with silicon gasket Press fit Style	Straight Type S.S. 304 Nozzle
 <p data-bbox="204 1809 469 1865">POLY CARBONATE 500 ML SQUARE SHAPE WITH SS CAP MODEL NO: SS-50C-BT</p>	SS50C-BT	500ml Square Shape	Polycarbonate Vacuum Style Without Ball	S.S. Cap with silicon gasket Press fit Style	Bent Type S.S. 304 Nozzle







Image	Model No	Capacity	Bottle Type Vacuum Style Without Ball	Cap	Sipper tube type
<b>Polysulfone Bottle for Rabbit / Guinea Pig</b>					
<p>POLYSULFONE 500 ML SQUARE SHAPE WITH SS CAP SS-50S-ST</p>	SS-50S-ST	500ml Square Shape	Polysulfone Vacuum Style Without Ball	S.S. Cap with silicon gasket Press fit Style	Straight Type S.S. 304 Nozzle
<p>POLYSULFONE 500 ML SQUARE SHAPE WITH SS CAP MODEL NO: SS-50S-BT</p>	SS50C-ST	500 ml Square Shape	Polycarbonate Vacuum Style Without Ball	S.S. Cap with silicon gasket Press fit Style	Straight Type S.S. 304 Nozzle

**Note:** For all above bottles of 650ml and 500ml capacity are kindly select Vacuum type or Ball Type Option while placing orders.

\*Autoclavable at 121°C temperature and pressure 15 LB/ Sq inch for 15 minutes.



## Orchid Scientific & Innovative India Pvt. Ltd.

-  B-59, M.I.D.C., Ambad, Nashik - 422010, India.
-  +91253-2387600, 2972525
-  office@orchidscientific.com, exports@orchidscientific.com
-  www.orchidscientific.com

