

Euroseries Rabbit Cage

EURO RABBIT CAGE SYSTEM

Euro Rabbit Cage System is the latest upgraded model with improved ergonomics, and reduced costs all while providing a safe, secure, and long-lasting housing solution for our partners and their research. This system features of large and comfortable spaces with the provision of removable partitions providing a facility for socializing for cage inhabitants.



Upgraded Model With Improved Animal Welfare, Ergonomics & Quality.



ANIMAL WELFARE

- Noryl perforated floor assures no sharp edges enhancing animal comfort.
- Smooth and rounded extrusions adjacent to the perforations allow more efficient drain into to excreta tray and also provide warm, quiet walk space reducing sore animal hocks.
- The rack interlinking system gives provision to increase floor space and socializing for cage inhabitants

ERGONOMICS AND QUALITY:

- Noryl polymer ensures repeated and longer sterilization cycles
- A fast-dry roof design results in quicker and easier washing procedures
- A push-and-close system allows the gate to be closed easily even while holding animals reducing stress in the animals.
- Plastic cage body resistant to alkaline and acid detergents
- Transparent removable side panels allow clear visibility of the entire cage

TECHNICAL DETAILS

Rack: The rack main frame is fabricated from AISI 304 quality stainless steel rectangular pipe of Size 30 x 20mm. Racks comprising of cage top & and side parts are made of S.S. 304 rods 2.5, 4.0, and 6.0 mm in diameter & and S.S. 304 sheets of 0.6, 1.2, 1.5 mm.

Door: The door of size 746mm x 385mm is fitted on the main frame of the trolley.

Floor Tray: Perforated floor tray of Noryl material of size 710 x 714 mm with 5/8"/16 m circular floor perforations for urine and feces collection.

Excreta Tray: 720 x720 mm Stainless Steel 304 for waste collection

Food Hopper: Fabricated from 23swg stainless steel 304 sheets of size 150x190x205mm and fits on the door.

Water Bottle: Made of autoclavable food grade Polycarbonate material of capacity 500ml with S.S. nozzle

Wheels: The entire battery with cages will move on 4 nos. of Polyurethane swivel 4" castor wheels two with brakes and two without brakes.

Side walls: Cage side walls are made of transparent polycarbonate material with openable by sliding. The middle partition in the case of 6 Cage system is also made of transparent polycarbonate material.

Rack Interlinking Clamps: One side of each rack has interlinking clamps for connecting or adding additional racks for floor space expansions, breeding, or socializing.

All the joints will be argon arc welded with a smooth finish.

Cage Free Height: 465mm

Cage Floor Space: 4200 sq. cm

Single Cage Inner dimensions: 725 x 715 x 530mm

ORDERING INFORMATION:

Product	Model	Outer Dimensions approx with wheel	Accessories Supplied
Euro Rabbit Cage System (6 cages)	RBINT-NOR-01 (6 cages)	1545 X 765 X 1810 mm	<ul style="list-style-type: none">• Trolley with sides, top and door fitted with castor wheels along with following accessories:• Noryl perforated trays: 06• S.S 304 waste collection trays: 06• Food Hopper: 06• Water bottle 500ml Polycarbonate: 06
Euro Rabbit Cage System (3 cages)	RBINT-NOR-02 (3 cages)	785 X 765 X 1810 mm	<ul style="list-style-type: none">• Trolley with sides, top and door fitted with castor wheels along with following accessories:• Noryl perforated trays: 03• S.S 304 waste collection trays: 03• Food Hopper: 03• Water bottle 500ml Polycarbonate: 03

Optional: Automatic water drinking nozzles with piping system and water tank.

**Needs to be specified in the order information*



Orchid Scientific & Innovative India Pvt. Ltd.

📍 B-59, M.I.D.C., Ambad, Nashik - 422010, India.

☎ +91253-2387600, 2972525

✉ office@orchidscientific.com, exports@orchidscientific.com

🌐 www.orchidscientific.com



www.orchidscientific.com