

Guinea Pig Pans: Model GPP-1

- Made of AISI non magnetic alkaline and chemical resistant 304 quality stainless steel material.
- Pans are made of S.S 304 grade heavy duty 16 SWG sheet.
- Pans are provided with 2 handles made from 4-5mm solid S.S. wire on each side for easy lifting
- Pan has provision for holding 2 water bottles and supplied with 2 number of P.P. 650ml bottles with S.S. nozzle.
- Each pan has holes of 20mm diameters on all sides
- All the joints will be argon arc welded with smooth and glossy finish.
- Overall inner size: 920mm x 600mm x 230mm (LXWXH)
- Floor Space: 5520 sq.cm



ORDERING INFORMATION:

Product	Model	Accessories Supplied
Guinea Pig Pan S.S.	GPP-1	• S.S. 304 Pan with 650ml bottle and bottle holder
OPTIONAL:	J Shape Food Hopper	

TROLLEY/ RACK FOR GUINEA PIG PAN :MODEL GPP-T-1

- Made from 25mm x25mm x 16swg S.S. 304 square pipe
- Suitable to hold 3 guinea pig pans with one pan in each shelf.
- Movable on 4 No. of Polyurethane swivel castor wheels 2 with brakes.
- The entire rack will be provided with stoppers on the backside of the angles to prevent cage falling from the back.
- All the joints will be welded with argon arc welding with glossy finish

ORDERING INFORMATION:

Product	Number of cages/shelf	Total cages/trolley
Rack for Guinea Pig Pan	GPP-T-1	3
Rack for Guinea Pig Pan	GPP-T-2	4



Orchid Scientific & Innovative India Pvt. Ltd.

- 📍 B-59, M.I.D.C., Ambad, Nashik - 422010, India.
- ☎ +91253-2387600, 2972525
- ✉ office@orchidscientific.com, exports@orchidscientific.com
- 🌐 www.orchidscientific.com

