

Animal Cages for Laboratory Small Rodents (Mice)

Orchid's Animal Cages are designed to facilitate animal welfare, meet research requirements and minimise the experimental variations. However, while selecting the cage, guidelines for laboratory animal facilities by concerned government organization regarding the minimum space requirements for housing of laboratory animals must be followed. As caging equipment is one of the most important factors in the physical and social environment of the experimental animal, Orchid has developed the largest range of cages meeting current norms and guidelines.

FEATURES:

- Wide range of animal house products for different animal species
- Durable and autoclavable (steam)
- Practicality of use
- Ergonomics

POLYPROPYLENE MICE CAGE (MODEL SMP 01:

Overall Size: 290 x 220 x 140mm (L X B X H)

Floor Space: 388cm²

CE Certified



BODY

- Autoclavable*
- Food Grade quality material Polypropylene material
- Uniform thickness of 4mm±0.5mm
- Excellent stackability for ease of handling and storage



**Designed To Facilitate
Animal Welfare, Meet
Research Requirements
And Minimise The
Experimental Variations**

S S TOP GRILL

- Made of S.S. 304 material with duly spot-welded and finished joints
- AISI 304 quality stainless steel material
- Provision for fixing water bottle
- Provision for Feed Pallets
- Fitted with two hooks of brass material for locking the grill to the top of the body of cage
- Fitted with tag holder
- Main frame is made of 4mm S.S. 304 rod and rest of grill of 2mm S.S. 304 rod
- Gap between S.S. rods is around 7mm
- Designed to ensure perfect fit and locking throughout the complete cage body
- Grill are designed such as to fit in second grill for ease of stacking with foldable partition for food and water bottle compartment

POLYCARBONATE MICE CAGE (SMC01):

Overall Size: 290 x 220 x 143mm (L X B X H)

Floor Space: 388cm²

CE Certified



BODY

- Autoclavable*
- Fully transparent body
- Food Grade quality material Polycarbonate material
- Uniform thickness of 4mm±0.5mm
- Excellent stackability for ease of handling and storage

S S TOP GRILL

- Made of S.S. 304 material with duly spot-welded finished joints and no sharp edges.
- AISI 304 quality stainless steel material.
- Provision for fixing water bottle.
- Provision for Feed Pallets.
- Fitted with two hooks of brass material for locking the grill to the top of the body of cage.
- Fitted with tag holder.
- Main frame is made of 4mm S.S. 304 rod and rest of grill of 2mm S.S. 304 rod.
- The gap between S.S. rods is around 7mm.
- Designed to ensure perfect fit and locking throughout the complete cage body.
- Grill are designed such as to fit in second grill for ease of stacking with foldable partition for food and water bottle compartment.

SPECIFICATION

Animal weight in gram	Floor Space of cage cm ²		Number of Mice per cage as per norms (CPCSEA)	
	SMP 01	SMC 01	SMP 01	SMC 01
<10	388	388	10	10
upto 15	388	388	7	7
upto 25	388	388	5	5
>25	388	388	4	4

ORDERING INFORMATION

Model	Suitable for animal	Overall Dimensions mm (L X B X H)	Components Included
SMP 01	Mice	290 x 220 x 140	P.P. Tray with S.S. grill
SMC 01	Mice	290 x 220 x 140	P.P. Tray with S.S. grill

*Autoclavable at 121°C temperature and pressure 15 LB/ Sq inch for 15 minutes.



Orchid Scientific & Innovative India Pvt. Ltd.

📍 B-59, M.I.D.C., Ambad, Nashik - 422010, India.

☎ +91253-2387600, 2972525

✉ office@orchidscientific.com, exports@orchidscientific.com

🌐 www.orchidscientific.com

