

Laboratory Animal Enrichment Devices

FOR ENHANCING ANIMAL WELFARE IN ANY SPECIES

Enrichment products designed to accommodate the behavior traits of individual animal species can stimulate lab animals' physical and mental activity. A research environment can be a stressful and unnatural setting for laboratory animals. Enrichment device enhances the animal's physiological and psychological well-being, prevents stereotypical behaviour, promotes the expression of species-specific behavior, and improves animal welfare. Therefore, Environmental Enrichment is an essential component of a successful research program. Environmental enrichment is an easy and efficient way of enhancing animal welfare in any species. As the importance of improving animal welfare has grown in significance, so has the range of Enrichment products.



Designed To Accommodate The Behaviour Traits Of Individual Animal Species Can Stimulate Lab Animals' Physical And Mental Activity

ENRICHMENT DEVICES FOR MOUSE AND RAT:

- Biohut for rat
- Biohut for mouse
- Biotunnel for rat
- Biotunnel for mouse
- Crawl-Ball
- Rat Tunnel
- Mouse Tunnel
- Mouse Fort

ENRICHMENT DEVICES FOR RABBIT & GUINEA PIG:

- Rabbit Huts
- Guinea-Pig Huts
- Jingle Ball
- Dumbbells

BIOHUT FOR RAT:

Biohut for Rat is an economical, rat enrichment product made of 100% pure virgin pulp paper. It gives a rat the opportunity to hide and also gives the rat a choice! Some rats choose to shred the huts for burrowing or nesting, while others shelter inside them.



BIOHUT FOR MOUSE:

- Economical and disposable enrichment
- Acrylamide Free
- Made of non-recycled, 100% virgin wood pulp
- 2 openings for increased visibility
- Separates easily to make 2 half huts
- Provides sheltering and nesting opportunities



BIOTUNNEL FOR RAT:

- Disposable
- Very economical sheltering, nesting and gnawing enrichment
- Made of 100% recycled pulp paper
- Helps reduce boredom and aggressive behavior
- Increases physical activity
- Not autoclavable contains adhesives



BIOTUNNEL FOR MOUSE:

- Disposable
- Very economical sheltering, nesting and gnawing enrichment
- Made of 100% recycled pulp paper
- Helps reduce boredom and aggressive behavior
- Increases physical activity
- Not autoclavable contains adhesives



CRAWL-BALL:

- Made of high quality polycarbonate
- Provides a comfortable hiding place
- Transparent material allows for ease of visibility
- Autoclavable & Cage washer safe



RAT TUNNEL:

- Made of high quality polycarbonate
- Offers rodents shelter and security
- Transparency makes routine cage checks easy and fast
- Helps reduce boredom and aggressive behavior
- Autoclavable and Cage washer safe
- To increase physical activity and usable floor space, suspend with the optional stainless steel hanger in the cage



MOUSE TUNNEL:

- Made of high-quality polycarbonate
- Offers rodent shelter and security
- Transparency makes routine cage checks easy and fast
- Helps reduce boredom and aggressive behavior
- Autoclavable and Cage washer safe
- To increase physical activity and usable floor space, suspend with the optional stainless steel hanger in the cage




MOUSEFORT:

- Made of high quality polycarbonate
- Offers rodents shelter and security
- Transparency makes routine cage checks easy and fast
- Low profile fits in most shoebox cages
- Provides shelter and a secure place to build a nest
- Autoclavable and Cage washer safe
- Angled sides prevent pushing against the water source




GUINEA-PIG HUTS:

- Made of polycarbonate
- Provides a comfortable hiding place
- Transparent material allows for ease of visibility
- Stacks easily for storage
- Cagewasher safe



RABBIT HUTS:

- Made of polycarbonate
- Provides a comfortable hiding place
- Transparent material allows for ease of visibility
- Stacks easily for storage
- Cage washer safe




JINGLE BALLS:

- Made of sturdy polyethylene
- Durable hollow ball, solid one piece construction
- 3 stainless steel "jingle" noise makers inside
- Slot is sealed to help prevent washer removal
- Helps to satisfy natural curiosity to explore and manipulate objects
- Cage washer safe



DUMBBELLS:

- Made of solid polypropylene
- Large ends knobs minimize loss through cage bars and floor slats
- Autoclavable & Cage washer safe



ORDERING INFORMATION:

Product Name	Packing	Material	Dimensions
Biohut for rat	80pcs/Box	Made of non-recycled virgin wood pulp	7 7/8" wide x 3 3/4" tall x 6" long
Biohut for mouse	500pcs/box	Made of non-recycled virgin wood pulp	3 1/4" wide x 2 1/2" tall x 5 1/4" long
Biotunnel for rat	100pcs/box	Made of non-recycled virgin wood pulp	6" long x 3" inside diameter x 1/12" wall
Biotunnel for mouse	200pcs/box	Made of non-recycled virgin wood pulp	3 7/8" long x 2" inside diameter; 1/12" wall
Crawls Balls	20pcs/box	Polycarbonate, Red colour	4" diameter with three 2 1/2" openings
Rat Tunnel	50pcs/box	Polycarbonate, Red color	6" long x 3" inside diameter; 1/12" wall
Mouse Tunnel	50pcs/box	Polycarbonate, Red color	3 7/8" long x 2" outside diameter; 1/12" wall
Mouse Fort	20pcs/box	Polycarbonate, Red colour	4.13" x 4.13" x 1.96" Four opening
Rabbit Huts	10pcs/box	Polycarbonate, Red colour	11" wide (ID) x 10" high x 12" long
Guinea-Pig Huts	10pcs/box	Polycarbonate, Red colour	6 1/2" wide (ID) x 4 1/2" high x 8" long
Jingle Balls	10pcs/box	Made of Polyethylene, Orange color	4 1/2" diameter
Dumbbells Small	10pcs/box	Made of P.P, Red/ Blue color	3 3/4" long with 2" ends
Dumbbells Medium	10pcs/box	Made of P.P, Red/ Blue color	4 3/4" long with 2 1/4" ends



Orchid Scientific & Innovative India Pvt. Ltd.

📍 B-59, M.I.D.C., Ambad, Nashik - 422010, India.

☎ +91253-2387600, 2972525

✉ office@orchidscientific.com, exports@orchidscientific.com

🌐 www.orchidscientific.com



www.orchidscientific.com