

Laboratory Animal Enrichment Devices

WOOD ENRICHMENT DEVICES


Orchid's Wood Enrichment Devices improves general well-being of Laboratory Animals & helps to maintain their dental health. These devices improve 3D area inside cage & are compatible for use with conventional cages and IVC systems. All devices are made from 100% Natural, Untreated, Specially Grown Certified Hardwood.

Enrichment Devices are manufactured from trees are specially grown in pollution-free land of Himachal Pradesh (India) without use of synthetic chemical fertilizers and pesticides, with manufacturing & plantation site being approx. 1400 meters above sea level & about 200 square kilometers away from any industrial set up. Our specially grown enrichment hardwood does not have aromatic hydrocarbons like softwood and are free from sharp nails or toxic chemical adhesive used for fixing. Orchid's Wood Enrichment Devices are compatible for use with both IVC and conventional cages.

Features:

- Made from 100% Natural, Untreated, Specially Grown Certified Hardwood
- Fulfills Natural Behaviors like Biting/Gnawing, Territory Marking, etc.
- Optimizes Feed Intake by avoiding unusual gnawing habits on feed pallet
- Improves General Well-Being & Maintain Dental Health of Lab Animals
- Gives complete Natural Environment & Improves 3D area inside the cage
- Most Economical & Suitable Choice as per Species-Specific Behaviors
- Most Safe. No risk of dioxins/BPA like with paper/plastic enrichment
- Autoclavable, Non-Toxic, Kiln Dried, Bio-degradable
- Material Safety: Batch Specific COA for screening of contaminants (heavy metals, mycotoxins, pesticides, PCBs, phytoestrogen, microbial load, etc.) tested in an NABL Accredited Lab




ORDERING INFORMATION:

Sr.No	Item	Picture	Unit Details	Pack Size
1	XS-Brick (Mice)		Enrichment Chew Sticks 50x10x10 mm	500pc
2	S-Brick (Mice)		Enrichment Chew Sticks 50x20x20 mm	500pc
3	M-Brick (Rat/Mice)		Enrichment Chew Sticks 100x20x20 mm	500pc
4	L-Brick (Rabbit/G.Pig)		Enrichment Chew Sticks 100x45x45 mm	200pc



Designed To Accommodate The Behaviour Traits Of Individual Animal Species Can Stimulate Lab Animals' Physical And Mental Activity

ORDERING INFORMATION:

Sr.No	Item	Picture	Unit Details	Pack Size
5	Cube (Mice)		Structural Enrichment 60x60x60 mm, Ø35 mm	50pc
6	Crossway (Mice)		Structural Enrichment 150x10x45 mm, Ø35 mm	50pc
7	Cube (Rat/Hamster)		Structural Enrichment 90x90x90 mm, Ø60 mm	50pc
8	Crossway (Rat/Hamster)		Structural Enrichment 170x15x75 mm, Ø55 mm	50pc



Orchid Scientific & Innovative India Pvt. Ltd.

📍 B-59, M.I.D.C., Ambad, Nashik - 422010, India.

☎ +91253-2387600, 2972525

✉ office@orchidscientific.com, exports@orchidscientific.com

🌐 www.orchidscientific.com



www.orchidscientific.com