

Trolley With Inbuilt Polypropylene Rabbit Cages

MODEL: RBINT-PP-01

The system has six rabbit cages and each rabbit cage comprises of P.P. tray forming bottom tray, S.S. door, top & side mesh inbuilt fitted on the trolley.



POLYPROPYLENE TRAY

The cage tray is made of autoclavable food grade P.P. material and has provision for fixing food hopper. P.P. Tray has provision for keeping floor net made of S.S. material. Size of P.P. tray is 650 x 500 x 215mm.

RACK

The rack structure will be fabricated from AISI 304 quality stainless steel square pipe of Size 30mm x 30mm. Sides of racks comprising cage top & side parts are made of 2.5mm diameter S.S. 304 rod.

DOOR

The door is fitted on the main frame of trolley and its size is 400mm x 195mm.

FLOOR GRILL:

The floor grill is made from 12 swg S.S. wire of size 12mm x 26mm welded on a frame of 5mm.

EXCRETA TRAY

Fabricated from 22swg AISI 304 quality stainless steel sheet. Each tray is fitted with handle in the front to pull the tray from trolley.

FOOD HOPPER

Fabricated from 22swg stainless steel 304 sheets of size 130 x 100 x 145mm and fits on the polypropylene body.

WATER BOTTLE

Made of autoclavable food grade P.P. material of capacity 650ml with S.S. nozzle

WHEELS

The entire battery with cages will move on 4 nos. of Polyurethane swivel 4" castor wheels two with brakes and two without brakes.

All the joints will be argon arc welded with smooth and glossy finish.

OVERALL SIZE

1432 x 525 x 1880 (L x B x H) with wheels

ORDERING INFORMATION:

Product	Model	Outer dimensions approx with wheel (mm)
Trolley with Inbuilt Rabbit cage	RBINT-PP-01 (6 cages)	1432 X 525 X 1880
Trolley with Inbuilt Rabbit cage	RBINT-PP-02 (3 cages)	731 x 525 x 1880



Orchid Scientific & Innovative India Pvt. Ltd.

📍 B-59, M.I.D.C., Ambad, Nashik - 422010, India.

☎ +91253-2387600, 2972525

✉ office@orchidscientific.com, exports@orchidscientific.com

🌐 www.orchidscientific.com



www.orchidscientific.com