

Rabbit Cage (Standard Type) Model RBS-S-01



MOC

Made of AISI non magnetic alkaline and chemical resistant 304 quality stainless steel sheet and wire.

MAIN FRAME

The main frame of the cage made of heavy duty S.S. square pipe of 12.7mm x 12.7mm x 16 swg and covered up to a height of 155mm from bottom with SS sheet of 22swg on three sides except front. Wire mesh made of 2.5mm and 5mm S.S. rod welded inside the SS square pipe on all three side including top of cage.

DOOR

Cage door has a fixed label holder and arrangement for fixing water bottle and food hopper. The door will be made of wire mesh of 2.5mm and 5mm SS rod welded on a frame of 5mm SS rod and covered with 22swg SS sheet.

FOOD HOPPER

J-shaped food hopper made of SS304 22swg

WATER BOTTLE

Polypropylene water bottle of capacity 650ml fitted S.S. bent type nozzle

FLOOR GRILL

Bottom floor net of 12mm x 26mm is made of 2.5mm & 5mm SS wire welded on a frame of 5mm SS rod. Floor net is fixed and non-removable type.

TRAY

Excreta Tray made of 22swg Stainless steel sheet with turned and smoothed edges. The cage is provided with S.S. L shaped angle of 20mm x20mm x 16swg for Excreta Tray.

CAGE SIZE

600mm x 450mm x 450mm (L x W x H)

FLOOR SPACE

2444 sq. cm

TROLLEY/ RACK TO HOLD RABBIT CAGES (STANDARD TYPE)

Available in Two options:

Sliding type and Hanging type.

Suitable for holding Rabbit Cage (Standard Type) Model RBS-S-01.

Material:

AISI 304 quality stainless steel square pipe of Size 30mm x 30mm x 16swg & stainless steel L shaped angle of 14swg. Movable on 4 No. of Polyurethane swivel castor wheels 2 with brakes. The entire Battery/Rack will be provided with stoppers on the backside of the angles to prevent cage falling from the back. All the joints will be welded with Argon Arc welding with glossy finish.

Rack will have inter-shelf gap of 21 inches in between two shelves so that each cage when placed in rack will have an extra space.

ORDERING INFORMATION

Model	Number of shelves	Number of cages/shelf	Total cages/trolley
RBS-S-T-3	3	1	3
RBS-S-T-6	3	2	6
RBS-S-T-9	4	3	9

Note: Trolleys for Rabbit cages of other sizes are also available on request.

RABBIT CAGE (BREEDING TYPE) : MODEL RBB-S-01



The breeding type cage basically comprises of 2 rabbit cages of size 600mm x 450mm x 450mm separated by a removable partition. The partition will be made of wire mesh of 2.5mm & 5mm SS rod welded on a frame of 5mm SS rod and is provided with a handle for ease in removing and inserting.

MOC

Made of AISI non magnetic alkaline and chemical resistant 304 quality stainless steel sheet and wire.

MAIN FRAME

Main frame of the cage made of heavy duty S.S. square pipe of 12.7mm x 12.7mm x 16 swg and covered up to a height of 155mm from bottom with SS sheet of 22swg on three sides except front. Wire mesh made of 2.5mm and 5mm S.S. rod welded inside the SS square pipe on all three side including top of cage.

DOOR

Cage has two doors and each door has a fixed label holder and arrangement for fixing water bottle and food hopper. The door will be made of wire mesh of 2.5mm and 5mm SS rod welded on a frame of 5mm SS rod and covered with 22swg SS sheet.

FOOD HOPPER

J shaped food hopper made of SS304 22swg

WATER BOTTLE

Polypropylene water bottle of capacity 650ml fitted with cap having S.S. bent type nozzle

FLOOR GRILL

Each part is provided with one bottom floor net. Bottom floor net is of 12x26mm made of 2.5mm & 5mm SS wire welded on a frame of 5mm SS rod. Floor net is fixed and non-removable type.

TRAY

Excreta tray made of 22swg Stainless steel sheet with turned and smoothed edges. The cage is provided with S.S. L shaped angle of 20mmx20mmx16swg for Excreta Tray.

CAGE SIZE

1210mm x 450mm x 450mm (L x W x H)

ORDERING INFORMATION

Product Name	Model	Accessories Supplied
Rabbit Cage Standard Model	RBS-S-01	<ul style="list-style-type: none">• S.S. 304 cage with floor grill• Food hopper• Excreta tray• Water bottle
Rabbit Cage Breeding Type	RBB-S-01	<ul style="list-style-type: none">• S.S. 304 cage with floor grill• Food hopper• Excreta tray• Water bottle

TROLLEY / RACK TO HOLD RABBIT CAGES (BREEDING TYPE):

Available in two options:

Sliding type and Hanging type.

Suitable for holding Rabbit Cage (Breeding Type) Model RBB-S-01 of size 1210mm x 450mm x 450mm (L x W x H)

Material:

AISI 304 quality stainless steel square pipe of Size 30mm x 30mm x 16swg & stainless steel L shaped angle of 14swg. Movable on 4 No. of Polyurethane swivel castor wheels 2 with brakes.

The entire Battery/Rack will be provided with stoppers on the backside of the angles to prevent the cage falling from the back. All the joints will be welded with Argon Arc welding with glossy finish.

Rack will have inter-shelf gap of 21 inches in between two shelves so that each cage when placed in rack will have an extra space.

ORDERING INFORMATION

Model	Number of shelves	Number of cages/shelf	Total cages/trolley
RBB-S-T-3	3	1	3



Orchid Scientific & Innovative India Pvt. Ltd.

📍 B-59, M.I.D.C., Ambad, Nashik - 422010, India.

☎ +91253-2387600, 2972525

✉ office@orchidscientific.com, exports@orchidscientific.com

🌐 www.orchidscientific.com

