

TROLLEY WITH INBUILT POLYPROPYLENE GUINEA PIG CAGES

MODEL GPINT-PP-01

System has six or eight guinea pig cages and each guinea pig cage comprises of P.P. tray forming bottom tray, S.S. door, top & side mesh inbuilt fitted on the trolley.



POLYPROPYLENE TRAY

The Tray is made of autoclavable food grade P.P. material and has provision for fixing food hopper. Size of P.P. tray is 650 x 500 x 215mm.

RACK

The rack structure will be fabricated from AISI 304 quality stainless steel square pipe of Size 30mm x 30mm. Sides of racks comprising cage top & side parts are made of 2.5mm diameter S.S. 304 rod.

DOOR

The door is fitted on the main frame of trolley and its size is 400mm x 195mm.

FOOD HOPPER

Fabricated from S.S. 304 sheet of size 130 x 100 x 145mm and fits on the polypropylene body.

WATER BOTTLE

Made of autoclavable food grade P.P. material of capacity 650ml with S.S. nozzle.

WHEELS

The entire battery with cages will move on 4 nos. of Polyurethane swivel 4” castor wheels two with brakes and two without brakes.

All the joints will be argon arc welded with smooth and glossy finish.

OVERALL SIZE

1432 x 525 x 1600 (L x B x H) without wheels

ORDERING INFORMATION:

Product	Model	Accessories Supplied
Trolley with inbuilt 6 P.P. guinea pig cage	GPINT-PP-01	Trolley with sides, top and door fitted, G.P. Cage tray- 6, each cage supplied with food hopper and water bottles
Trolley with inbuilt 8 P.P. guinea pig cage	GPINT-PP-02	Trolley with sides, top and door fitted, G.P. Cage tray- 8, each cage supplied with food hopper and water bottles



Orchid Scientific & Innovative India Pvt. Ltd.

📍 B-59, M.I.D.C., Ambad, Nashik - 422010, India.

☎ +91253-2387600, 2972525

✉ office@orchidscientific.com, exports@orchidscientific.com

🌐 www.orchidscientific.com

