

Trolleys For Mice And Rat Cages

Racks / Trolleys for holding animal cages either mice or rat are available in 2 types- Hanging Type and Sliding Type.

In Hanging Type, the Cage is kept in the rack with the collar of cage in hanging position on the supports provided in the rack.

In case of Sliding Type, the cage is kept in the rack with its base sliding on the supports provided in the rack.



HANGING TYPE RACKS FOR MICE CAGES:

Suitable for holding Polypropylene or Poly-carbonate Cages of size 29 x 22 x 14 cm. Made of AISI 304 quality S.S. square pipe of Size 25 mm x 25 mm and S. S. 25 x 3mm angle for hanging the cage. Movable on 4 No. Of Polyurethane swivel 4" castor wheels out of which 2 are with brakes. All the joints will be welded with Argon Arc welding with glossy /mat finish.

Model	Suitable for Cage size	No of Cages/row	Total Shelves	Total Number of Cages
ACTH-M01	29X22X14cm	5	3	15
ACTH-M02	29X22X14cm	5	5	25

SLIDING TYPE RACKS FOR MICE CAGES:

Suitable for holding Polypropylene or Poly-carbonate Cages of size 29 x 22 x 14cm. Made of AISI 304 quality stainless steel square pipe of Size 25 mm x 25 mm and stainless steel angle. Movable on 4 No. of Polyurethane swivel 4" castor wheels 2 with brakes. The entire Rack is provided with stoppers on the backside of the angles to prevent the falling of cage from back. All the joints will be welded with Argon Arc welding with Glossy / Mat finish.

Model	Suitable for Cage size	No of Cages/row	Total Shelves	Total Number of Cages
ACTH-M01	29X22X14cm	5	3	15
ACTH-M02	29X22X14cm	5	5	25
ACTS-M02	29X22X14cm	6	6	36

HANGING TYPE RACKS FOR RAT CAGES:

Suitable for holding Polypropylene or Poly-carbonate Cages of size 43 x 29 x 20 cm. or 41 x 28 x 15 cm.* Made of AISI 304 Quality S. S. square pipe of Size 25 mm x 25mm and S. S. 25 x 3mm angle for hanging the cage. Movable on 4 No. of Polyurethane swivel 4" castor wheels out of which 2 are with brakes. All the joints will be welded with Argon Arc welding with glossy /mat finish.

Model	Suitable for Cage size	No of Cages/row	Total Shelves	Total Number of Cages
ACTH-R01	43X29X20/41X28X15cm	5	3	15
ACTH-R02	43X29X20/41X28X15cm	5	5	25
ACTH-R03	65x50x20cm	4	4	16

SLIDING TYPE RACKS FOR RAT CAGES:

Suitable for holding Polypropylene or Poly-carbonate Cages of size 43 x 29 x 20 cm. or 41 x 28 x 15 cm.* Made of AISI 304 quality stainless steel square pipe of Size 25mm x 25 mm and stainless steel angle for sliding the cage. Movable on 4 No. of Polyurethane swivel 4" castor wheels 2 with brakes. The entire Rack is provided with stoppers on the backside of the angles to prevent the falling of cage from back. All the joints will be welded with Argon Arc welding with Glossy / Mat finish.



Model	Suitable for Cage size	No of Cages/row	Total Shelves	Total Number of Cages
ACTS-R01	43X29X20/41X28X15cm	5	3	15
ACTS -R02	43X29X20/41X28X15cm	5	5	25
ACTS -R03	43X29X20cm	6	6	36
ACTS-R04	65 X 50 X 20CM	4	4	16

All Sliding Type of Racks are supplied in parts ready to assemble and are supplied with required fasteners and tools.

All Hanging type of Racks are in intact form and are not supplied in parts, so while ordering ease of transportation and packing must be taken into consideration.



Orchid Scientific & Innovative India Pvt. Ltd.

📍 B-59, M.I.D.C., Ambad, Nashik - 422010, India.

☎ +91253-2387600, 2972525

✉ office@orchidscientific.com, exports@orchidscientific.com

🌐 www.orchidscientific.com

