

GUINEA PIG PANS : MODEL GPP-1

- Made of AISI non magnetic alkaline and chemical resistant 304 quality stainless steel material.
- Pans are made of S.S 304 grade heavy duty 16 SWG sheet.
- Pans are provided with 2 handles made from 4-5mm solid S.S. wire on each side for easy lifting.
- Pan has provision for holding 2 water bottles and supplied with 2 number of P.P. 650ml P.P. bottles with S.S. nozzle.
- Each pan has holes of 20mm diameters on all sides.
- All the joints will be argon arc welded with smooth and glossy finish.
- Overall inner size: 920mm x 600mm x 230mm (LXWXH).
- Floor Space : 5520 cm²



ORDERING INFORMATION

PRODUCT	MODEL	ACCESSORIES SUPPLIED
GUINEA PIG PAN S. S.	GPP-1	S.S. 304 Pan with 650ml bottle and bottle holder
OPTIONAL : J Shape Food Hopper		

TROLLEY/ RACK FOR GUINEA PIG PAN : Model GPP-T-1

- Made from 25mm x 25mm x 16swg S.S. 304 square pipe.
- Suitable to hold 3 guinea pig pans with one pan in each shelf.
- Movable on 4 No. of Polyurethane swivel castor wheels 2 with brakes.
- The entire rack will be provided with stoppers on the backside of the angles to prevent cage falling from the back.
- All the joints will be welded with argon arc welding with glossy finish.

ORDERING INFORMATION

PRODUCT	MODEL	NUMBER OF PANS
RACK FOR GUINEA PIG PAN	GPP-T-1	3
RACK FOR GUINEA PIG PAN	GPP-T-2	4



ORCHID SCIENTIFIC & INNOVATIVE INDIA PVT. LTD.

B 59, M. I. D. C., AMBAD, NASHIK- 422010

TEL: +91253-2387600, 2382525, 2381515 FAX: +91253-2382525

W : www.orchidscientifics.com E : orchidscientific@gmail.com, info@orchidscientifics.com

Orchid's continuing products improvements program makes specifications subject to change without notice