

QUAIL CAGE

The Japanese quail (*Coturnix japonica*) has been a popular animal model in numerous fields of research. Orchid has designed special cages for housing of Quail suiting to laboratory and research application.



Features

- Made from complete AISI S.S. 304 quality material
- Overall Dimensions : 900 X 600 X 300mm (L X W X H)
- Frame is made from 12.7 X 12.7 X 1.5mm S.S. 304 square pipe
- Food Hopper : Supplied with two J shaped food hoppers of 22swg.
- Water drinking: Supplied with specially designed water system of dimensions 300 X 50 X 50mm made from 1.5mm thick S.S. 304 sheet material.
- Door: Door is fitted with label holder and has dimensions of 343 X 275mm.
- Top, front and side grills are made from 3mm & 5mm S.S. 304 bars and have gap of 10 X 43mm.
- Both sides are fitted with handles for ease of holding during the shifting of cages.
- Floor grill is made from 2.5mm and 3mm S.S. 304 bar with gap of 10x10mm.
- Tray dimensions: 928 X 597 X 75mm, Made from S.S. 304 sheet of 0.8mm thickness.
Tray has to be placed in the trolley below the mounted cage.
- Floor Space : 5031 sq.cm

Note: Orchid's continuing product development makes specifications subject to change without prior notification.



ORCHID SCIENTIFIC & INNOVATIVE INDIA PVT. LTD.

B 59, M. I. D. C., AMBAD, NASHIK- 422010

TEL: +91253-2387600, 2382525, 2381515 FAX: +91253-2381515

W : www.orchidscientific.com

E : orchidscientific@gmail.com, info@orchidscientific.com

Orchid's continuing products improvements program makes specifications subject to change without notice

QUAIL CAGE



Trolley for Quail Cage :

- Sliding type trolley suitable to hold 3 cages of size 900 X 600 X 300mm along with trays.
- Made from S.S.304 square pipe of 25 X 25 X 1.5mm and angle of 25 X 25 X 2mm.
- Overall dimensions: 995 X 650 X 1600mm.
- Trolley is supplied with four castor wheels out of which two are with brake and two are without brakes.

ORDERING INFORMATION :

Product Name	Model	Accessories Supplied
Quail Cage	QC-01	Cage with floor grill-1, Water drinking system-1, Food hopper-02 Excreta tray-1
Trolley for Quail cages	QC-T-01	Trolley with 4 castor wheels

Note: Orchid's continuing product development makes specifications subject to change without prior notification.



ORCHID SCIENTIFIC & INNOVATIVE INDIA PVT. LTD.

B 59, M. I. D. C., AMBAD, NASHIK- 422010

TEL: +91253-2387600, 2382525, 2381515 FAX: +91253-2381515

W : www.orchidscientific.com

E : orchidscientific@gmail.com, info@orchidscientific.com

Orchid's continuing products improvements program makes specifications subject to change without notice