

TROLLEY WITH INBUILT POLYPROPYLENE GUINEA PIG CAGE : MODEL GPINT-PP-01

System has six or eight guinea pig cages and each guinea pig cage comprises of P. P. tray forming bottom tray, S. S. door, top & side mesh inbuilt fitted on the trolley.

- **Polypropylene Tray** : Tray is made of autoclavable food grade P. P. material and has provision for fixing food hopper. Size of P. P. tray is 650 x 500 x 215mm.
- **Rack** : The rack structure will be fabricated from AISI 304 quality stainless steel square pipe of Size 30mm x 30mm. Sides of racks comprising cage top & side parts are made of 2.5mm diameter S.S. 304 rod.
- **Door** : The door is fitted on the main frame of trolley and its size is 400mm x 195mm.

- **Food Hopper** : Fabricated from S.S. 304 sheet of size 130x100x145mm and fits on the polypropylene body.
- **Water Bottle** : Made of autoclavable food grade P.P. material of capacity 650ml with S.S. nozzle.
- **Wheels** : The entire battery with cages will move on 4 nos. of Polyurethane swivel 4" castor wheels two with brakes and two without brakes.
- All the joints will be argon arc welded with smooth and glossy finish.
- **Overall Size** : 1432 x 525 x 1600mm (L x B x H) without wheels.

ORDERING INFORMATION

PRODUCT	MODEL	ACCESSORIES SUPPLIED
Trolley With Inbuilt 6 P. P. Guinea Pig Cage	GPINT-PP-01	Trolley with sides, top and door fitted, G. P. Cage tray- 6, each cage supplied with food hopper and water bottles
Trolley With Inbuilt 8 P. P. Guinea Pig Cage	GPINT-PP-02	Trolley with sides, top and door fitted, G. P. Cage tray- 8, each cage supplied with food hopper and water bottles



ORCHID SCIENTIFIC & INNOVATIVE INDIA PVT. LTD.

B 59, M. I. D. C., AMBAD, NASHIK- 422010

TEL: +91253-2387600, 2382525, 2381515 FAX: +91253-2382525

W : www.orchidscientifics.com E : orchidscientific@gmail.com, info@orchidscientifics.com

Orchid's continuing products improvements program makes specifications subject to change without notice